

# Goat Anti-Human IgG+IgM+IgA H&L ab102416

## 1 References

### 概述

产品名称	山羊抗人IgG+IgM+IgA H&L
宿主	Goat
靶标种属	Human
特异性	By immunoelectrophoresis and ELISA this antibody reacts with human IgG, IgA and IgM. No antibody was detected against non-immunoglobulin serum proteins. This antibody may cross react with immunoglobulins from other species.
经测试应用	适用于: ICC, IHC-P, ELISA, WB
免疫原	Other Immunogen Type corresponding to IgG + IgM + IgA.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 6.8 Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
纯化说明	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads.
克隆	多克隆
同种型	IgG

### 应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab102416于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC		Use at an assay dependent dilution.

应用	Ab评论	说明
IHC-P		1/200 - 1/2000.
ELISA		1/100 - 1/30000. ELISA - Coating 1/100 - 1/500; ELISA - Primary 1/1000 - 1/30000
WB		1/1000 - 1/30000. WB - Colorimetric 1/1000 - 1/10000; WB - Chemiluminescent 1/1000 - 1/30000

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors