

Goat Anti-Monkey IgG H&L (HRP) ab112767

9 References

概述

产品名称	山羊Anti-Monkey IgG H&L (HRP)
宿主	Goat
靶标种属	Monkey
特异性	By immunoelectrophoresis and ELISA ab112767 reacts specifically with Monkey IgG and with light chains common to other Monkey immunoglobulins. Species detected are Rhesus, Baboon, Carb Green, Cyno, Chimp and African Green. Reaction to human was also detected. No antibody was detected against non-immunoglobulin serum proteins. May cross react with IgG from other species.
经测试应用	适用于: WB, IHC, ICC, ELISA
免疫原	Other Immunogen Type corresponding to IgG.
偶联物	HRP

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.80 Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative, 99% PBS
纯度	Immunogen affinity purified
纯化说明	ab112767 was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to horseradish peroxidase (HRP).
共轭说明	Molar enzyme/antibody protein ratio is 4:1.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab112767于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/5000 - 1/50000. Predicted molecular weight: 36 kDa.
IHC		1/200 - 1/5000.
ICC		1/200 - 1/5000.
ELISA		1/10000 - 1/100000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors