

### Goat Anti-Mouse IgE (SureLight® APC) ab99904

#### 概述

产品名称	山羊抗小鼠IgE (SureLight® APC)
宿主	Goat
靶标种属	Mouse
经测试应用	适用于: Flow Cyt
偶联物	SureLight® APC. Ex: 652nm, Em: 657nm

#### 性能

形式	Lyophilized: Lyophilized material is stable for up to one year. Reconstitute with 1.0 ml of ddH <sub>2</sub> O. After product has been reconstituted, product should be stored at 2-8°C in the dark and be used within 3 months. If further dilution of the conjugate is required, use diluted material within one week.
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 8.20 Constituents: 0.87% Sodium chloride, 0.121% Tris, 1.71% Sucrose, 0.1% BSA
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

#### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用 ab99904 于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度; 实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use a concentration of 3 µg/ml.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors