

Goat Anti-Mouse IgG Fc (6nm Gold) preadsorbed ab105285

2 References

概述		
产品名称	山羊抗小鼠IgG Fc (6nm Gold)预吸附二抗	
宿主	Goat	
靶标种属	Mouse	
经测试应用	适用于: Electron Microscopy, WB, ELISA	
最小交叉反应	Cow, Horse, Human	
偶联物	Gold 6nm	more details

性能		
形式	Liquid	
存放说明	Shipped at 4°C. Store at +4°C.	
存储溶液	pH: 8.50 Preservative: 0.05% Sodium azide Constituents: 1.5% BSA, 0.164% Sodium phosphate, 0.381% Sodium borate, 0.87% Sodium chloride	
纯度	Immunogen affinity purified	
纯化说明	purified by affinity chromatography	
克隆	多克隆	
同种型	IgG	

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab105285于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Electron Microscopy		1/20 - 1/40.

应用	Ab评论	说明
WB		Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors