

Goat Anti-Mouse IgG Fc (FITC) ab97264

5 References 1 图像

概述	
产品名称	山羊抗小鼠IgG Fc (FITC)
宿主	Goat
靶标种属	Mouse
特异性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Mouse IgG-Fc. Cross reactivity with other immunoglobulins and light chains is less than 0.1%.
经测试应用	适用于: IHC-P, ICC/IF
偶联物	FITC. Ex: 493nm, Em: 528nm

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, PBS
纯度	Immunogen affinity purified
纯化说明	Antiserum was solid phase adsorbed to ensure class specificity. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to FITC.
共轭说明	F/P ratio is 4.4
克隆	多克隆
同种型	IgG

应用


The Abpromise guarantee **Abpromise™**承诺保证使用ab97264于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/50 - 1/500.

应用	Ab评论	说明
ICC/IF		1/50 - 1/500.

图片



A431 cell nucleus (DAPI signal) showing three gamma-H2AX foci structures.

Primary antibody H2A.X (phospho S139) antibody [3F2], [ab22551](#), 100ug.

Secondary antibody Mouse IgG-Fc (FITC), ab97264, 1mg.

Immunocytochemistry/ Immunofluorescence - Goat

Anti-Mouse IgG Fc (FITC) (ab97264)

This image is courtesy of Jorge E. Gonzalez (CPHR, La Habana, Cuba), Joan F. Barquinero (UAB, Barcelona, Spain) and Jessica Martinez (UAB, Universitat Autònoma de Barcelona, Spain).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors