

Product datasheet

Goat Anti-Mouse IgG H&L (Alkaline Phosphatase) preadsorbed ab7069

★★★★★ 1 Abreviews 4 References 1 图像


概述	
产品名称	山羊抗小鼠IgG H&L (Alkaline Phosphatase)预吸附二抗
宿主	Goat
靶标种属	Mouse
经测试应用	适用于: ICC/IF, ELISA, IHC-P, IHC-Fr, Immunomicroscopy, Dot blot, WB
最小交叉反应	Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Human, Rabbit, Rat, Sheep more details
免疫原	Full length native mouse IgG (purified).
偶联物	Alkaline Phosphatase

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium azide Constituents: 0.00136% Zinc chloride, 0.0095% Magnesium chloride, 0.79% Tris HCl, 50% Glycerol (glycerin, glycerine), 0.87% Sodium chloride, 1% BSA
纯度	Affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using mouse IgG coupled to agarose beads.
共轭说明	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
克隆	多克隆

应用	
The Abpromise guarantee	Abpromise™ 承诺保证使用ab7069于以下的经测试应用
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。	

应用	Ab评论	说明
ICC/IF		1/200 - 1/1000.
ELISA		1/7000.
IHC-P	★★★★★ (1)	Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.
Immunomicroscopy		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution.
WB		1/1000 - 1/2000.

图片



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Goat Anti-Mouse IgG H&L (Alkaline Phosphatase) preadsorbed (ab7069)
This image is courtesy of an anonymous Abreview.

ab13918 staining RANK in human skin tissue sections by Immunohistochemistry (IHC-P - paraformaldehyde-fixed, paraffin-embedded sections). Tissue was fixed with formaldehyde and blocked with 1% BSA/5% normal goat serum/0.1% saponin for 30 minutes at 22°C. Samples were incubated with primary antibody (1/50 in blocking buffer) for 9 hours at 4°C. An alkaline phosphatase-conjugated goat anti-mouse IgG H&L (ab7069) (1/300) was used as the secondary antibody.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors