

Goat Anti-Mouse IgG H&L (DyLight® 550) preadsorbed ab96880

[10 References](#) [1 图像](#)

概述

产品名称	山羊抗小鼠IgG H&L (DyLight® 550)预吸附二抗
宿主	Goat
靶标种属	Mouse
特异性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Mouse IgG and with light chains common to other Mouse immunoglobulins. No antibody was detected against non immunoglobulin serum proteins. Reduced cross-reactivity to bovine, chicken, goat, horse, human, pig, rabbit and rat IgG was detected.
经测试应用	适用于: IHC-P, ICC/IF, Flow Cyt
最小交叉反应	Chicken, Cow, Goat, Horse, Human, Pig, Rabbit, Rat more details
偶联物	DyLight® 550. Ex: 562nm, Em: 576nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, PBS
纯度	Immunogen affinity purified
纯化说明	Antiserum was cross adsorbed using bovine, chicken, horse, human, pig, rabbit and rat immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 550.
克隆	多克隆
同种型	IgG
常规说明	DyLight® 550 replaces DyLight® 549.

应用


The Abpromise guarantee

Abpromise™ 承诺保证使用 ab96880 于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/50 - 1/500. DyLight® 550 replaces DyLight® 549. The DyLight® 550 excitation is 562nm and the emission is 576nm. DyLight® 550 provides the same orange-to-red fluorescence and photostability as the DyLight® 549 dye. Please refer to the vial in order to check whether the product you received is DyLight® 549 or DyLight® 550. This information however will have no impact on your
ICC/IF		1/50 - 1/500.
Flow Cyt		1/50 - 1/200.

图片



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors