



Product datasheet

Goat Anti-Mouse IgG H&L (PE) preadsorbed ab99911

概述

产品名称	山羊抗小鼠IgG H&L (PE)预吸附二抗	
宿主	Goat	
靶标种属	Mouse	
经测试应用	适用于: Flow Cyt	
最小交叉反应	Cow, Horse, Human, Pig, Rabbit	<a href="#">more details</a>
偶联物	PE. Ex: 488nm, Em: 575nm	

性能

形式	Lyophilized:Please reconstitute in 1mL of distilled water, vortex gently and wait for twenty minutes. The effective concentration after re suspension should be 0.1 mg/mL.
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at +4°C. Store In the Dark.
存储溶液	pH: 7.40 Preservative: 0.05% Sodium azide Constituents: 0.1% BSA, 0.87% Sodium chloride, 1.64% Sodium phosphate, 6.84% Sucrose
纯度	Protein G purified
共轭说明	R-phycoerythrin Excitation max.:565>498 nm Emission max: 578 nm
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee      Abpromise™承诺保证使用ab99911于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use a concentration of 3 µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors