

# Goat Anti-Mouse IgG1 heavy chain (Biotin) preadsorbed ab98691

## 9 References

### 概述

产品名称	山羊抗小鼠IgG1 heavy chain (Biotin)预吸附二抗	
宿主	Goat	
靶标种属	Mouse	
经测试应用	适用于: IHC-P, ICC/IF, ELISA, WB	
最小交叉反应	Human, Rat	<a href="#">more details</a>
偶联物	Biotin	

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA
纯度	Immunogen affinity purified
纯化说明	ab98691 was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Biotin.
克隆	多克隆
同种型	IgG

### 应用

**The Abpromise guarantee** **Abpromise™** 承诺保证使用ab98691于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/250 - 1/2500.

应用	Ab评论	说明
ICC/IF		1/100 - 1/500.
ELISA		1/1000 - 1/10000. as primary
WB		1/10000 - 1/200000. by chemiluminescent detection

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors