

Goat Anti-Mouse IgG+IgM H&L (Biotin) preadsorbed ab47844


3 References 1 图像

概述	
产品名称	山羊抗小鼠IgG+IgM H&L (Biotin)预吸附二抗
宿主	Goat
靶标种属	Mouse
特异性	This antibody reacts with Mouse IgG and IgM heavy and light chains. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, anti-Biotin, Mouse IgG and Mouse IgM.
经测试应用	适用于: ELISA, WB, ICC/IF, IHC-P
最小交叉反应	Cow, Horse, Human
免疫原	Full length IgG and IgM molecules (Mouse)
偶联物	Biotin
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
纯度	Immunogen affinity purified
纯化说明	Prepared from polyspecific antiserum by immunoaffinity chromatography using antigens coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
共轭说明	10-20 BAC (Biotinamidocaproate N-Hydroxysuccinimide Ester) molecules per IgG molecule
克隆	多克隆
同种型	IgG
应用	

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		1/20000 - 1/100000.
WB		1/2000 - 1/10000.
ICC/IF		1/1000 - 1/5000.
IHC-P		Use at an assay dependent dilution.

## 图片



Ab47844 was used at dilution 1/1000 with the primary antibody **ab24575** in ICC/IF. See the review on **ab24575**.

Immunocytochemistry/ Immunofluorescence - Goat  
Anti-Mouse IgG+IgM H&L (Biotin) preadsorbed  
(ab47844)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors