

Goat Anti-Mouse IgG+IgM H&L (HRP) preadsorbed ab47827

26 References 1 图像

概述	
产品名称	山羊抗小鼠IgG+IgM H&L (HRP)预吸附二抗
宿主	Goat
靶标种属	Mouse
特异性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, anti-Peroxidase, Mouse IgG and Mouse IgM.
经测试应用	适用于: ICC/IF, ELISA, WB, IHC-P, IHC-Fr, Electron Microscopy
最小交叉反应	Cow, Horse, Human
免疫原	Full length native proteins (purified) (Mouse)
偶联物	HRP

[more details](#)

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.50 Preservative: 0.01% Gentamicin sulphate Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
纯度	IgG and protease free Immunogen affinity purified
纯化说明	Prepared from polyspecific antiserum by immunoaffinity chromatography using antigens coupled to agarose beads followed by solid phase adsorption (s) to remove any unwanted reactivities.
克隆	多克隆
同种型	IgG


应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab47827于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		1/500 - 1/2500.
ELISA		1/5000 - 1/40000.
WB		1/2000 - 1/10000.
IHC-P		1/500 - 1/2500.
IHC-Fr		1/500 - 1/2500.
Electron Microscopy		Use at an assay dependent dilution.

图片



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Goat Anti-Mouse IgG+IgM H&L (HRP) preadsorbed (ab47827)

Ab47827 was used at dilution 1/500 with the primary antibody **ab2869** in IHC-P. See the review on **ab2869**.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors