

### Goat Anti-Rabbit IgG H&L (DyLight® 488) ab96883

★★★★★ [2 Abreviews](#) [26 References](#) [1 图像](#)

#### 概述

|       |  |
|-------|--|
| 产品名称  | 山羊抗兔IgG H&L (DyLight® 488)   |
| 宿主    | Goat   |
| 靶标种属  | Rabbit   |
| 特异性   | By immunoelectrophoresis and ELISA this antibody reacts specifically with Rabbit IgG and with light chains common to other Rabbit immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. |
| 经测试应用 | <b>适用于:</b> ICC/IF, IHC-P, Flow Cyt  |
| 免疫原   | Other Immunogen Type corresponding to Rabbit IgG.  |
| 偶联物   | DyLight® 488. Ex: 493nm, Em: 518nm   |

#### 性能

|      |  |
|------|--|
| 形式   | Liquid   |
| 存放说明 | Shipped at 4°C. Store at +4°C.   |
| 存储溶液 | pH: 6.8<br>Preservative: 0.09% Sodium azide<br>Constituents: 0.2% BSA, PBS   |
| 纯度   | Immunogen affinity purified  |
| 纯化说明 | This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 488. |
| 克隆   | 多克隆  |
| 同种型  | IgG  |


#### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab96883于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用       | Ab评论      | 说明            |
|----------|-----------|---------------|
| ICC/IF   |           | 1/50 - 1/500. |
| IHC-P    | ★★★★★ (1) | 1/50 - 1/500. |
| Flow Cyt |           | 1/50 - 1/200. |

## 图片



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors