

Goat Anti-Rabbit IgG H&L (HRP) preadsorbed ab7090

★★★★★ [2 Abreviews](#) [552 References](#) [2 图像](#)

概述	
产品名称	山羊抗兔IgG H&L (HRP)预吸附二抗
宿主	Goat
靶标种属	Rabbit
经测试应用	适用于: ICC/IF, IHC-P, WB, ELISA, Immunomicroscopy, Dot blot, IHC-Fr
最小交叉反应	Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Human, Mouse, Rat, Sheep <a href="#">more details</a>
免疫原	Full length native Rabbit IgG (purified).
偶联物	HRP

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.50 Preservative: 0.01% Gentamicin sulphate Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
纯度	Affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose beads.
共轭说明	Horseradish Peroxidase (HRP)
克隆	多克隆
同种型	IgG

应用	
The Abpromise guarantee	<a href="#">Abpromise™</a> 承诺保证使用ab7090于以下的经测试应用
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。	

应用	Ab评论	说明
ICC/IF		1/1000 - 1/5000.
IHC-P	★★★★★ (1)	Use at an assay dependent concentration.
WB	★★★★★ (1)	Use at an assay dependent concentration.
ELISA		1/20000 - 1/80000.
Immunomicroscopy		Use at an assay dependent concentration.
Dot blot		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.

图片

vehicle

OE

IPE

C

DSS


Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Goat Anti-Rabbit IgG H&L (HRP) preadsorbed (ab7090)

Lean QY et al. PLoS One (2015) Reproduced under the Creative Commons license  
<http://creativecommons.org/licenses/by/4.0/>

Immunohistochemical analysis of healthy and colitis mouse colon sections (untreated and treated with enoxaparin) labeling claudin-4 with **ab15104** at 1/200. ab7090, anti-rabbit immunoglobulin G conjugated to horseradish peroxidase (HRP) at 1/300 was used as the secondary antibody. **ab64238** at 1/50 was used to develop the histological signal.

Antigen retrieval was performed by incubating the sections for 10 minutes at 97°C in 1 mM EDTA buffer, pH 8.0 or 10 mM citrate buffer, pH 6.0.

Control, C; untreated colitis, DSS; oral enoxaparin, OE; intraperitoneal injection of enoxaparin, IPE.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Goat Anti-Rabbit IgG H&L (HRP) preadsorbed (ab7090)

This image is courtesy of an anonymous Abreview.

**ab9986** staining osteoprotegerin in human skin tissue sections by Immunohistochemistry (IHC-P - paraformaldehyde-fixed, paraffin-embedded sections). Tissue was fixed with formaldehyde and blocked with 1% BSA and 5% normal goat serum in PBS for 30 minutes at 22°C; Samples were incubated with primary antibody (1/20) for 9 hours at 4°C. An HRP-conjugated goat anti-rabbit IgG H&L (ab7090) (1/300) was used as the secondary antibody.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors