

# Goat anti-Rabbit IgG (SureLight® APC) ab72567

## 2 References

### 概述

产品名称	山羊抗兔IgG (SureLight® APC)
宿主	Goat
靶标种属	Rabbit
经测试应用	适用于: Flow Cyt, FRET
最小交叉反应	Cow, Mouse
免疫原	Rabbit IgG
偶联物	SureLight® APC. Ex: 652nm, Em: 657nm

[more details](#)

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.40 Preservative: 0.013% Sodium azide Constituents: 0.1% BSA, 0.87% Sodium chloride, 1.64% Sodium phosphate, 1.71% Sucrose
纯度	Size exclusion
纯化说明	Size exclusion chromatography as well as some proprietary methods for purifying the final materials were used.
共轭说明	Excitation max. $\lambda$ 652 nm. Emission max. $\lambda$ 657.5 nm. F/P 1:1.
克隆	多克隆
同种型	IgG

### 应用

**The Abpromise guarantee**

**Abpromise™** 承诺保证使用ab72567于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration.
FRET		Use at an assay dependent concentration.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors