

Goat Anti-Rat IgG H&L (DyLight® 488) preadsorbed ab98420

5 References 2 图像

概述		
产品名称	山羊抗大鼠IgG H&L (DyLight® 488)预吸附二抗	
宿主	Goat	
靶标种属	Rat	
经测试应用	适用于: Flow Cyt, ICC/IF, IHC-P	
最小交叉反应	Chicken, Cow, Goat, Horse, Human, Mouse, Rabbit, Sheep	
偶联物	DyLight® 488. Ex: 493nm, Em: 518nm	
		more details

性能		
形式	Liquid	
存放说明	Shipped at 4°C. Store at +4°C.	
存储溶液	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA	
纯度	Immunogen affinity purified	
纯化说明	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 488.	
克隆	多克隆	
同种型	IgG	

应用


The Abpromise guarantee **Abpromise™**承诺保证使用ab98420于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		1/50 - 1/200.

应用	Ab评论	说明
ICC/IF		1/50 - 1/500.
IHC-P		1/50 - 1/500.


图片



Spectrum - Goat Anti-Rat IgG H&L (DyLight® 488)
preadsorbed (ab98420)

Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.



Immunocytochemistry/ Immunofluorescence - Goat
Anti-Rat IgG H&L (DyLight® 488) preadsorbed
(ab98420)

ICC/IF image of **ab6161** stained HeLa cells. The cells were 100% methanol fixed (5 min) and then incubated in 1%BSA / 10% normal goat serum / 0.3M glycine in 0.1% PBS-Tween for 1h to permeabilise the cells and block non-specific protein-protein interactions. The cells were then incubated with the antibody (**ab6161**, 1µg/ml) overnight at +4°C. The secondary antibody (green) was DyLight® 488 goat anti-rat- H&L, pre-adsorbed (ab98420) used at a 1/250 dilution for 1h. Alexa Fluor® 594 WGA was used to label plasma membranes (red) at a 1/200 dilution for 1h. DAPI was used to stain the cell nuclei (blue) at a concentration of 1.43µM.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors