

Goat Anti-Rhesus monkey IgG H&L (PE) ab125852

★★★★★ 1 Abreviews

概述

| | |
|-------|---|
| 产品名称 | 山羊抗恒河猴IgG H&L (PE) |
| 宿主 | Goat |
| 靶标种属 | Rhesus monkey |
| 经测试应用 | 适用于: IHC-P, Dot blot, WB, ELISA, ICC/IF |
| 免疫原 | Full length native protein (purified) corresponding to IgG. Pooled antisera from goats hyperimmunized with rhesus IgG |
| 偶联物 | PE. Ex: 488nm, Em: 575nm |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C. Store In the Dark. |
| 存储溶液 | pH: 7.3 |
| 纯度 | Affinity purified |
| 纯化说明 | Affinity chromatography on rhesus IgG covalently linked to agarose. |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab125852于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|----------|------|--|
| IHC-P | | Use at an assay dependent concentration. |
| Dot blot | | Use at an assay dependent concentration. |
| WB | | Use at an assay dependent concentration. |

| 应用 | Ab评论 | 说明 |
|--------|------|--|
| ELISA | | Use at an assay dependent concentration. |
| ICC/IF | | Use at an assay dependent concentration. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors