



Product datasheet

Goat Anti-Syrian Hamster IgG H&L (Biotin) preadsorbed ab7145

1 References

| | | |
|--------|-------------------------------------------------------------------------|------------------------------|
| 概述 | | |
| 产品名称 | 山羊抗叙利亚仓鼠IgG H&L (Biotin)预吸附二抗 | |
| 宿主 | Goat | |
| 靶标种属 | Syrian hamster | |
| 经测试应用 | 适用于: ICC/IF, Dot blot, ELISA, IHC-P, WB, Immunomicroscopy, IHC-Fr | |
| 最小交叉反应 | Chicken, Cow, Goat, Guinea pig, Horse, Human, Mouse, Rabbit, Rat, Sheep | more details |
| 免疫原 | hamster IgG whole molecule | |
| 偶联物 | Biotin | |

| | | |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
| 性能 | | |
| 形式 | Liquid | |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. | |
| 存储溶液 | pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.88% Sodium chloride, 1% BSA, 0.42% Tripotassium orthophosphate | |
| 纯度 | Affinity purified | |
| 纯化说明 | This product was prepared from monospecific antiserum by immunoaffinity chromatography using hamster IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. | |
| 共轭说明 | Biotinamidocaproate N-Hydroxysuccinimide Ester (BAC) Biotin/Protein Ratio: 10-20 BAC molecules per Rabbit IgG molecule | |
| 克隆 | 多克隆 | |

应用

The Abpromise guarantee Abpromise™承诺保证使用ab7145于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|------------------|------|----|
| ICC/IF | | |
| Dot blot | | |
| ELISA | | |
| IHC-P | | |
| WB | | |
| Immunomicroscopy | | |
| IHC-Fr | | |

应用说明 immunoblotting, ELISA, immunohistochemistry, immunomicroscopy as well as other antibody based assays using streptavidin or avidin conjugates requiring lot-to-lot consistency

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors