

# Goat Anti-Syrian hamster IgG H&L (HRP) preadsorbed ab102327

### 概述

|        |  |
|--------|--|
| 产品名称   | 山羊抗叙利亚仓鼠IgG H&L (HRP)预吸附二抗   |
| 宿主     | Goat   |
| 靶标种属   | Syrian hamster   |
| 特异性    | By immunoelectrophoresis and ELISA this antibody reacts specifically with syrian hamster IgG and with light chains common to other syrian hamster immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. Less than 1% cross reactivity to bovine, chicken, goat, horse, human, mouse, rabbit, rat and sheep IgG was detected. This antibody may cross react with IgG from other species. |
| 经测试应用  | 适用于: IHC-P, ELISA, WB, ICC   |
| 最小交叉反应 | Chicken, Cow, Goat, Horse, Human, Mouse, Rabbit, Rat, Sheep <a href="#">more details</a>   |
| 偶联物    | HRP  |

### 性能

|      |  |
|------|--|
| 形式   | Liquid   |
| 存放说明 | Shipped at 4°C. Store at +4°C.   |
| 存储溶液 | pH: 6.8<br>Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative   |
| 纯度   | Immunogen affinity purified  |
| 纯化说明 | Antiserum was cross adsorbed using bovine, chicken, goat, horse, human, mouse, rabbit, rat and sheep immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to HRP. |
| 克隆   | 多克隆  |
| 同种型  | IgG  |

### 应用

The Abpromise guarantee

[Abpromise™](#) 承诺保证使用ab102327于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用    | Ab评论 | 说明  |
|-------|------|---|
| IHC-P |      | 1/200 - 1/5000.   |
| ELISA |      | 1/10000 - 1/100000. (primary).  |
| WB    |      | 1/5000 - 1/30000. (colorimetric), 1/10000 - 1/50000 (chemiluminescent). |
| ICC   |      | Use at an assay dependent dilution.                                     |

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors