



Product datasheet

Anti-Human IgG4 pFc' [HP6023] ab99814

概述

产品名称	抗人IgG4 pFc' [HP6023]
宿主	Mouse
靶标种属	Human
特异性	ab99814 reacts with the pFc' portion of the heavy chain of Human IgG4.
经测试应用	<b>适用于:</b> ELISA, WB, Dot blot, IHC-Fr, ICC/IF
免疫原	Full length protein corresponding to IgG4.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 8.20 Constituent: 100% Borate buffered saline  No preservatives or amine-containing buffer salts added.
纯度	Protein A purified
克隆	单克隆
克隆编号	HP6023
同种型	IgG3
轻链类型	kappa

应用

The Abpromise guarantee      **Abpromise™**承诺保证使用ab99814于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.

应用	Ab评论	说明
WB		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution. (for Dot and slot immunoblotting)
IHC-Fr		Use at an assay dependent dilution.
ICC/IF		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors