

Mouse Anti-Human IgG4 pFc' [HP6023] (HRP) ab99817

13 References 1 图像

概述

产品名称	小鼠抗人IgG4 pFc' [HP6023] (HRP)
宿主	Mouse
靶标种属	Human
经测试应用	适用于: ICC, ELISA, Dot blot, WB, IHC-Fr
免疫原	Full length protein corresponding to Human IgG4. Human IgG4 myeloma protein
偶联物	HRP

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.40 Constituents: 50% Glycerol, 50% PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	HP6023
同种型	IgG3
轻链类型	kappa

应用


The Abpromise guarantee      **Abpromise™**承诺保证使用ab99817于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC		Use at an assay dependent dilution.
ELISA		1/1000 - 1/2000.

应用	Ab评论	说明
Dot blot		Use at an assay dependent dilution. for Dot and slot immunoblotting.
WB		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.

## 图片



ELISA plate was coated with purified human IgG1, IgG2, IgG3, and IgG4. Immunoglobulins were detected with serially diluted ab99817.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors