

Mouse monoclonal [4E3] Anti-Human IgG1 hinge heavy chain (Alkaline Phosphatase) ab99773

2 References

概述	
产品名称	小鼠单克隆抗体[4E3]抗人IgG1 hinge heavy chain (Alkaline Phosphatase)
宿主	Mouse
靶标种属	Human
特异性	ab99773 reacts with the hinge region of the heavy chain of Human IgG1 as demonstrated by ELISA. May also react with IgG from other species.
经测试应用	适用于: ELISA, Flow Cyt
偶联物	Alkaline Phosphatase

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 8.00 Preservative: 0.1% Sodium azide Constituents: 50% Glycerol (glycerin, glycerine), 0.0095% Magnesium chloride, 0.605% Tris
纯度	Protein G purified
克隆	单克隆
克隆编号	4E3
同种型	IgG1

应用

The Abpromise guarantee      **Abpromise™**承诺保证使用ab99773于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		1/500 - 1/1000.

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors