

## Product datasheet

# Mouse monoclonal [4E3] Anti-Human IgG1 hinge heavy chain (HRP) ab99774

★★★★★ [1 Abreviews](#) [17 References](#)

### 概述

|       |   |
|-------|---|
| 产品名称  | 小鼠单克隆抗体[4E3]抗人IgG1 hinge heavy chain (HRP)  |
| 宿主    | Mouse                                       |
| 靶标种属  | Human                                       |
| 经测试应用 | <b>适用于:</b> ELISA, ICC                      |
| 免疫原   | Other Immunogen Type corresponding to IgG1. |
| 偶联物   | HRP   |

### 性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Store at +4°C.                  |
| 存储溶液 | pH: 7.40<br>Constituents: 50% Glycerol, 50% PBS |
| 纯度   | Protein G purified                              |
| 克隆   | 单克隆   |
| 克隆编号 | 4E3   |
| 同种型  | IgG1  |

### 应用

**The Abpromise guarantee** **Abpromise™** 承诺保证使用ab99774于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用    | Ab评论      | 说明                                  |
|-------|-----------|-------------------------------------|
| ELISA | ★★★★★ (1) | 1/1000 - 1/2000.                    |
| ICC   |           | Use at an assay dependent dilution. |

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors