

# Mouse monoclonal [B3102E8] Anti-Human IgE Fc ab99804

## 1 图像

### 概述

|       |  |
|-------|--|
| 产品名称  | 小鼠单克隆抗体[B3102E8]抗人IgE Fc   |
| 宿主    | Mouse  |
| 靶标种属  | Human  |
| 特异性   | ab99804 reacts with the Fc portion of the heavy chain of Human IgE |
| 经测试应用 | 适用于: ELISA, WB, Dot blot, IHC-Fr, ICC                              |

### 性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.   |
| 存储溶液 | pH: 8.20<br>Constituent: 100% Borate buffered saline  |
|      | No preservatives or amine-containing buffer salts added.  |
| 纯度   | Affinity purified   |
| 克隆   | 单克隆   |
| 克隆编号 | B3102E8   |
| 同种型  | IgG1  |
| 轻链类型 | kappa   |
| 常规说明 | This product has switched from ascites to tissue culture supernatant on 8 <sup>th</sup> June 2020. Lot numbers higher than GR3258714 are from tissue culture supernatant. Please note that the dilutions may need to be adjusted accordingly. If you have any questions, please do not hesitate to contact our scientific support team. |


### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab99804于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用       | Ab评论 | 说明                                  |
|----------|------|-------------------------------------|
| ELISA    |      | Use at an assay dependent dilution. |
| WB       |      | Use at an assay dependent dilution. |
| Dot blot |      | Use at an assay dependent dilution. |
| IHC-Fr   |      | Use at an assay dependent dilution. |
| ICC      |      | Use at an assay dependent dilution. |

## 图片



ELISA plate was coated with Human IgE. Immunoglobulin was detected with serially diluted Mouse monoclonal [B3102E8] Anti-Human IgE Fc (ab99804) conjugated to Biotin followed by Streptavidin-HRP and ABTS.

This product has switched from ascites to tissue culture supernatant on 8<sup>th</sup> June 2020. Lot numbers higher than GR3258714 are from tissue culture supernatant.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors