

Product datasheet

Mouse monoclonal [B3102E8] Anti-Human IgE Fc (HRP) ab99806

★★★★★ [1 Abreviews](#) [9 References](#)

概述

产品名称	小鼠单克隆抗体[B3102E8]抗人IgE Fc (HRP)
宿主	Mouse
靶标种属	Human
经测试应用	适用于: ELISA, WB, Dot blot, IHC-Fr, ICC
偶联物	HRP

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.40 Constituents: 50% Glycerol, 50% PBS
纯度	Proprietary Purification
纯化说明	Purified from TCS. The unlabeled antibody is manufactured using a combination of precipitation and/or chromatography techniques yielding a product of = 90% purity by HPLC and/or SDS-PAGE.
克隆	单克隆
克隆编号	B3102E8
同种型	IgG1
轻链类型	kappa
常规说明	This product was changed from ascites to tissue culture supernatant on 13 July 2019. Please note that the dilutions may need to be adjusted accordingly. If you have any questions, please do not hesitate to contact our scientific support team.

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab99806于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA	★★★★★ (1)	Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.
Dot blot		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
ICC		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors