

Mouse monoclonal [B3506B4] Anti-Human IgA1 Fc (PE) ab99797

概述

产品名称	小鼠单克隆抗体[B3506B4]抗人IgA1 Fc (PE)
宿主	Mouse
靶标种属	Human
特异性	ab99797 reacts with the Fc portion of the heavy chain of IgA1.
经测试应用	适用于: ELISA, Dot blot, WB, IHC-Fr, ICC/IF
偶联物	PE. Ex: 488nm, Em: 575nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	pH: 7.3 Preservative: 0.1% Sodium azide Constituents: Sucrose, PBS
纯度	Protein G purified
克隆	单克隆
克隆编号	B3506B4
同种型	IgG1
轻链类型	kappa
常规说明	Source: Ascites fluid

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab99797于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution. for Dot and slot immunoblotting
WB		Use at an assay dependent dilution. Predicted molecular weight: 38 kDa.
IHC-Fr		Use at an assay dependent dilution.
ICC/IF		Use 0.2µg for 10 ⁶ cells.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors