

# Mouse monoclonal [hlgG1] Anti-Human IgG H&L preadsorbed ab436

## 1 References

### 概述

产品名称	小鼠单克隆抗体[hlgG1]抗人IgG H&L预吸附二抗
宿主	Mouse
靶标种属	Human
特异性	hlgG1 specifically recognizes human IgG. hlgG - 100% hlgG1 - 147% hlgG2 - 150% hlgG3 - 94% hlgG4 - 2.1% hlgM - 0.6% hlgA - 0.3% Goat IgG - <0.01% Rabbit IgG - <0.01% Horse IgG - <0.01% Pig IgG - <0.01% Cow IgG - <0.01% Sheep IgG - <0.01% Affinity Constant: $2 \times 10^{10}$ L/mol human IgG.
经测试应用	适用于: ELISA, RIA
最小交叉反应	Cow, Goat, Horse, Pig, Rabbit, Sheep

[more details](#)

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 6.00 Preservative: 0.095% Sodium azide Constituents: 0.335% PBS, 1.0878% Sodium citrate, 0.9% Sodium chloride
纯度	Protein A purified
克隆	单克隆
克隆编号	hlgG1
同种型	IgG2a
轻链类型	unknown

### 应用

The Abpromise guarantee

Abpromise™承诺保证使用ab436于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		
RIA		

#### 应用说明

ELISA: Use at an assay dependant dilution.

RIA: Use at an assay dependant dilution.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors