

Mouse monoclonal [JDC-10] Anti-Human IgG Fc (PE) ab99761

2 References

概述	
产品名称	小鼠单克隆抗体[JDC-10]抗人IgG Fc (PE)
宿主	Mouse
靶标种属	Human
特异性	Reacts with the Fc portion of the heavy chain of all subclasses of human IgG.
经测试应用	适用于: ICC/IF, ELISA, WB, Dot blot, IHC-Fr
偶联物	PE. Ex: 488nm, Em: 575nm

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	pH: 7.3 Preservative: 0.1% Sodium azide Constituents: Sucrose, PBS
纯度	Affinity purified
克隆	单克隆
克隆编号	JDC-10
同种型	IgG1
轻链类型	kappa

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab99761于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		Use 0.2µg for 10 ⁶ cells.
ELISA		Use at an assay dependent dilution.

应用	Ab评论	说明
WB		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution. + slot-immunoblotting: Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors