

# Mouse monoclonal [MARE-1] Anti-Rat IgE heavy chain (Biotin) ab11666

## 概述

产品名称	小鼠单克隆抗体[MARE-1]抗大鼠IgE heavy chain (Biotin)
宿主	Mouse
靶标种属	Rat
特异性	This antibody is specific for the epsilon chain of rat IgE and does not cross-react with other classes of rat immunoglobulin.
经测试应用	<b>适用于:</b> WB, ELISA
免疫原	Rat IR162, IR1016, IR2 and IR410 IgE myeloma proteins.
偶联物	Biotin

## 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium azide Constituents: PBS, 50% Glycerol
纯度	Immunogen affinity purified
纯化说明	Purified IgG prepared by affinity chromatography on IR2 rat IgE myeloma.
克隆	单克隆
克隆编号	MARE-1
骨髓瘤	Sp2/0
同种型	IgG1
常规说明	This product should be stored undiluted. Storage in frost free freezers is not recommended. Should this product contain a precipitate we recommend microcentrifugation before use.

## 应用

The Abpromise guarantee

**Abpromise™** 承诺保证使用ab11666于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use at an assay dependent dilution. Predicted molecular weight: 51.4 kDa.
ELISA		Use a concentration of 0.5 µg/ml.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors