



Product datasheet

Mouse monoclonal [UHB] Anti-Human IgM mu chain ab99741

概述	
产品名称	小鼠单克隆抗体[UHB]抗人IgM mu chain
宿主	Mouse
靶标种属	Human
特异性	ab99741 reacts with the heavy chain of Human IgM as demonstrated by ELISA; may also react with IgM from other species.
经测试应用	适用于: ELISA, Flow Cyt
免疫原	Other Immunogen Type corresponding to IgM.

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 8.20 Constituent: 100% Borate buffered saline
纯度	No preservatives or amine-containing buffer salts added. Protein A purified
克隆	单克隆
克隆编号	UHB
同种型	IgG3

应用	
The Abpromise guarantee	Abpromise™承诺保证使用ab99741于以下的经测试应用
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。	

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors