

Mouse monoclonal [UHB] Anti-Human IgM mu chain (Alkaline Phosphatase) ab99743

概述

产品名称	小鼠单克隆抗体[UHB]抗人IgM mu chain (Alkaline Phosphatase)
宿主	Mouse
靶标种属	Human
特异性	ab99743 reacts with the heavy chain of Human IgM as demonstrated by ELISA; may also react with IgM from other species.
经测试应用	适用于: Dot blot, WB, IHC-Fr, ELISA, ICC
免疫原	Other Immunogen Type corresponding to IgM.
偶联物	Alkaline Phosphatase

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 8.00 Preservative: 0.1% Sodium azide Constituents: 50% Glycerol, 0.0095% Magnesium chloride, 0.605% Tris
纯度	Protein A purified
克隆	单克隆
克隆编号	UHB
同种型	IgG3

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab99743于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Dot blot		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.
ELISA		1/500 - 1/1000.
ICC		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors