

# Mouse monoclonal [UHB] Anti-Human IgM mu chain (PE) ab99746

### 概述

产品名称	小鼠单克隆抗体[UHB]抗人IgM mu chain (PE)
宿主	Mouse
靶标种属	Human
特异性	ab99746 reacts with the heavy chain of Human IgM as demonstrated by ELISA. It may also react with IgM from other species.
经测试应用	<b>适用于:</b> ICC, Dot blot, ICC/IF, IHC-Fr, ELISA, WB
偶联物	PE. Ex: 488nm, Em: 575nm

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	pH: 7.3 Preservative: 0.1% Sodium azide Constituents: Sucrose, PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	UHB
同种型	IgG3

### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab99746于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution.

应用	Ab评论	说明
ICC/IF		Use 0.1 µg for 10 <sup>6</sup> cells. for Indirect immunofluorescent staining of IgM+ human B lymphocytes.
IHC-Fr		Use at an assay dependent dilution.
ELISA		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors