

Rabbit Anti-Armenian hamster IgG H&L (HRP) ab5745

★★★★★ 1 Abreviews 10 References 1 图像

概述

产品名称	兔抗亚美尼亚仓鼠IgG H&L (HRP)
宿主	Rabbit
靶标种属	Armenian hamster
特异性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti-Rabbit Serum, Armenian Hamster IgG and Armenian Hamster Serum.
经测试应用	适用于: Dot blot, WB, ELISA, Electron Microscopy, IHC-P, ICC
免疫原	Armenian Hamster IgG whole molecule.
偶联物	HRP

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.50 Preservative: 0.01% Gentamicin sulphate Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
纯化说明	BSA is IgG and protease free This product was prepared from monospecific antiserum by immunoaffinity chromatography using Armenian Hamster IgG coupled to agarose.
共轭说明	Horseradish Peroxidase (HRP).
克隆	多克隆
同种型	IgG


应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab5745于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Dot blot		Use at an assay dependent concentration.
WB		1/1000 - 1/5000.
ELISA		1/25000. This product has been assayed against 1.0 ug of Armenian Hamster IgG in a standard capture ELISA using ABTS (2,2'-azino-bis-[3-ethylbenthiiazoline-6-sulfonic acid]) as a substrate for 30 minutes at room temperature.
Electron Microscopy		Use at an assay dependent concentration.
IHC-P	★★★★★ (1)	1/1000 - 1/5000.
ICC		1/1000 - 1/5000.

图片



ab5745 tested against purified Armenian Hamster IgG. Each well was coated in duplicate with 1.0 µg of Armenian Hamster IgG. The starting dilution of antibody was 5µg/ml and the X-axis represents the log10 of a 3-fold dilution. This titration is a 4-parameter curve fit where the IC50 is defined as the titer of the antibody. Assay performed using 3% fish gelatin as blocking buffer and TMB substrate

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors