

## Product datasheet

# Rabbit F(ab')<sub>2</sub> Anti-Mouse IgG - H&L (Biotin), pre-adsorbed ab98775

### 概述

产品名称	兔F(ab') <sub>2</sub> 抗小鼠IgG - H&L (Biotin), 预吸附二抗
宿主	Rabbit
靶标种属	Mouse
经测试应用	适用于: IHC-P, ICC/IF, ELISA, WB
最小交叉反应	Human, Rat
偶联物	Biotin

[more details](#)

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA
纯度	Immunogen affinity purified
纯化说明	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. F(ab') <sub>2</sub> fragment were generated using a pepsin digestion. Fc fragments and whole IgG molecules have been removed. Fragments were conjugated to Biotin.
克隆	多克隆
同种型	IgG

### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab98775于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/500 - 1/5000.
ICC/IF		1/200 - 1/1000.
ELISA		1/20000 - 1/400000. (Primary).
WB		1/20000 - 1/400000. (Chemiluminescent).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors