

# Rabbit Anti-Horse IgG H&L (HRP) ab6921

## 9 References

### 概述

产品名称	兔抗马IgG H&L (HRP)
宿主	Rabbit
靶标种属	Horse
经测试应用	适用于: ICC/IF, IHC-Fr, Immunomicroscopy, WB, IHC-P, ELISA, Dot blot
免疫原	Full length native horse IgG (purified).
偶联物	HRP

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.01% Gentamicin sulphate Constituents: 0.878% Sodium chloride, 1% BSA, 0.424% Potassium phosphate
纯度	Affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using horse IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities
共轭说明	Horseradish Peroxidase (HRP)
克隆	多克隆
同种型	IgG

### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab6921于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		1/500 - 1/2500.

应用	Ab评论	说明
IHC-Fr		Use at an assay dependent dilution.
Immunomicroscopy		Use at an assay dependent dilution.
WB		1/5000 - 1/20000.
IHC-P		Use at an assay dependent dilution.
ELISA		1/50000 - 1/200000.
Dot blot		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors