

Rabbit Anti-Horse IgG+IgM+IgA H&L (HRP) ab8533

概述

产品名称	兔抗马IgG+IgM+IgA H&L (HRP)
宿主	Rabbit
靶标种属	Horse
特异性	Antibodies were raised in rabbit to Horse IgG, IgM and IgA
经测试应用	适用于: ELISA, Dot blot
偶联物	HRP

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.50 Constituents: 0.164% Sodium phosphate, 50% Glycerol (glycerin, glycerine), 0.58% Sodium chloride
纯度	Immunogen affinity purified
纯化说明	Specific antibodies were absorbed by incubation with Sepharose-bound Horse immunoglobulins. Specific antibodies were eluted by acidic buffer at pH 2.5 followed by neutralisation and dialysis. After repeated binding with immobilized Horse immunoglobulins a minimum of 65% protein bound.
共轭说明	The purified polyclonal was conjugated to horseradish peroxidase according to the periodate method followed by gel-filtration and ion-exchange chromatography to clear unbound conjugate.
克隆	多克隆

应用

The Abpromise guarantee      **Abpromise™**承诺保证使用ab8533于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		

应用	Ab评论	说明
Dot blot		

#### 应用说明

Recommended dilution as secondary antibody in ELISA is 1:1000.

Recommended dilution as secondary antibody in dot-blotting is 1:100.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors