

### Rabbit Anti-Human IgG H&L (HRP) preadsorbed ab98576

#### 概述

产品名称	兔抗人IgG H&L (HRP)预吸附二抗
宿主	Rabbit
靶标种属	Human
经测试应用	适用于: IHC-P, ELISA, WB, ICC
最小交叉反应	Mouse, Rat
偶联物	HRP

[more details](#)

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.8 Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative
纯度	Immunogen affinity purified
纯化说明	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to HRP.
克隆	多克隆
同种型	IgG

#### 应用

##### The Abpromise guarantee

**Abpromise™** 承诺保证使用ab98576于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/200 - 1/5000.
ELISA		1/10000 - 1/100000.

应用	Ab评论	说明
WB		1/5000 - 1/50000.
ICC		1/200 - 1/5000.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors