

Rabbit Anti-Human IgG4 H&L ab86251

1 References

概述	
产品名称	兔抗人IgG4 H&L
宿主	Rabbit
靶标种属	Human
经测试应用	适用于: WB, ELISA, RIA
免疫原	Full length native Immunoglobulin G4 purified from Human plasma
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: 49.98% PBS, 50% Glycerol (glycerin, glycerine)
纯度	Protein G purified
克隆	多克隆
同种型	IgG
应用	

The Abpromise guarantee      **Abpromise™**承诺保证使用ab86251于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use at an assay dependent dilution. Predicted molecular weight: 36 kDa.
ELISA		Use at an assay dependent dilution.
RIA		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors