

# Rabbit Anti-Human IgM H&L (FITC) ab8506

## 1 References

### 概述

产品名称	兔抗人IgM H&L (FITC)
宿主	Rabbit
靶标种属	Human
经测试应用	适用于: ICC/IF
偶联物	FITC. Ex: 493nm, Em: 528nm

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.50 Constituents: 0.164% Sodium phosphate, 50% Glycerol (glycerin, glycerine), 0.58% Sodium chloride
纯度	Immunogen affinity purified
纯化说明	From the crude polyclonal the crossreactive antibodies were extracted by incubation with Sepharose-bound human IgA and IgG. Specific antibodies were absorbed by incubation with Sepharose-bound human IgM. Specific antibodies were eluted by acidic buffer at pH 2.5 followed by neutralisation and dialysis. After repeated binding with immobilized human IgM a minimum of 65% protein bound.
共轭说明	The purified polyclonal was conjugated to fluorescein isothiocyanate isomer I followed by gel-filtration and ion-exchange chromatography to clear unbound conjugate.
克隆	多克隆
常规说明	F/P ratio is 4-6 mol Fluorescein per 1 mol of goat IgM. There is low background fluorescence in normal tissue in contrast to cases of autoimmune disease where IgM is deposited.

### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab8506于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		

**应用说明**

For immunohistochemical staining of human tissue sections a 1:40 dilution is recommended.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

**Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors