

Rabbit Anti-Mouse IgG H&L (Alkaline Phosphatase) ab6729

25 References

概述

产品名称	兔抗小鼠IgG H&L (Alkaline Phosphatase)
宿主	Rabbit
靶标种属	Mouse
经测试应用	适用于: Dot blot, ELISA, IHC-P, IHC-Fr, Immunomicroscopy, ICC/IF, WB
免疫原	Mouse IgG whole molecule
偶联物	Alkaline Phosphatase

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Diluted stock solutions should be used immediately and not be frozen. Do Not Freeze. Please see notes section. Store undiluted.
存储溶液	Preservative: 0.01% Sodium azide Constituents: 0.87% Sodium chloride, 1% BSA, 0.79% Tris HCl, 0.01% Magnesium chloride, 0.001% Zinc chloride, 50% Glycerol
纯度	Affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Mouse IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
共轭说明	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
克隆	多克隆
常规说明	Freezing alkaline phosphatase conjugates will result in a substantial loss of enzymatic activity.

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab6729于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Dot blot		Use at an assay dependent concentration.
ELISA		Use at an assay dependent dilution.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
Immunomicroscopy		Use at an assay dependent concentration.
ICC/IF		1/500 - 1/2000.
WB		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors