

Rabbit Anti-Mouse IgM mu chain ab9175

1 References

概述

产品名称	兔抗小鼠IgM mu chain
宿主	Rabbit
靶标种属	Mouse
特异性	The affinity purified anti IgM mu antibodies only react with IgM by immunodiffusion and Immunoelectrophoretic techniques. A 1 mg/ml solution of the antibody is reacted against normal mouse serum and IgG1, IgG2a, IgG2b, IgG3, IgM and IgA. By ELISA, there is minimal reactivity to IgG.
经测试应用	适用于: ELISA, Immunoelectrophoresis, WB
免疫原	Rabbits were immunized with highly purified normal mouse IgM in Freund's adjuvant.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.20
纯度	Immunogen affinity purified
纯化说明	The antiserum was solid phase absorbed to make it IgM mu chain specific and the anti mu chain antibodies were isolated by solid phase immunoaffinity chromatography.
克隆	多克隆

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab9175于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
Immunoelectrophoresis		Use at an assay dependent concentration.

应用	Ab评论	说明
WB		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors