

Rabbit Anti-Sheep IgM H&L (Alkaline Phosphatase) ab112761

1 References

概述

产品名称	兔抗羊IgM H&L (Alkaline Phosphatase)
宿主	Rabbit
靶标种属	Sheep
特异性	By immunoelectrophoresis and ELISA ab112761 reacts specifically with sheep IgM. May cross react with IgM from other species.
经测试应用	适用于: WB, ELISA, IHC-P, ICC/IF
免疫原	Other Immunogen Type corresponding to IgM.
偶联物	Alkaline Phosphatase

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, 0.58% Sodium chloride, 0.01% Magnesium chloride, 0.001% Zinc chloride, 1.19% HEPES
纯度	Immunogen affinity purified
纯化说明	Antiserum was solid phase adsorbed to ensure class specificity. ab112761 was isolated by affinity chromatography using antigen coupled to agarose beads.
共轭说明	Molar enzyme: antibody protein ratio is 1:1.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab112761于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/1000 - 1/10000. for chemiluminescent detection
ELISA		1/1000 - 1/10000. as a primary
IHC-P		1/100 - 1/1000.
ICC/IF		1/100 - 1/1000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors