

Rat monoclonal [23G3] Anti-Mouse IgE epsilon chain (HRP)  
ab99574

8 References

概述	
产品名称	大鼠单克隆抗体[23G3]抗小鼠IgE epsilon chain (HRP)
宿主	Rat
靶标种属	Mouse
特异性	ab99574 reacts with the epsilon heavy chain of BALB/c Mouse IgE.
经测试应用	适用于: ELISA, ICC/IF
偶联物	HRP
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.40 Constituents: 50% Glycerol, 50% PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	23G3
同种型	IgG1
轻链类型	kappa
应用	

The Abpromise guarantee      **Abpromise™**承诺保证使用ab99574于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		1/1000 - 1/2000.

应用	Ab评论	说明
ICC/IF		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors