

Rat monoclonal [JC5-1] Anti-Mouse lambda light chain (FITC)
ab99623

概述

产品名称	大鼠单克隆抗体[JC5-1]抗小鼠lambda light chain (FITC)
宿主	Rat
靶标种属	Mouse
特异性	ab99623 reacts with lambda 1 and lambda 2 light chains. It does not react with other immunoglobulin isotypes or kappa light chains.
经测试应用	适用于: ELISA, Flow Cyt, ICC/IF
偶联物	FITC. Ex: 493nm, Em: 528nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	pH: 7.4 Preservative: 0.1% Sodium azide Constituent: PBS
纯度	Affinity purified
克隆	单克隆
克隆编号	JC5-1
同种型	IgG2b
轻链类型	kappa

应用

The Abpromise guarantee Abpromise™承诺保证使用ab99623于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.

应用	Ab评论	说明
Flow Cyt		Use 1µg for 10 ⁶ cells. ab136125 - Rat monoclonal IgG2b, is suitable for use as an isotype control with this antibody.
ICC/IF		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors