

Rat monoclonal [LO-ME-3] Anti-Mouse IgE epsilon chain (Biotin) preadsorbed ab11580

7 References

概述

产品名称	大鼠单克隆抗体[LO-ME-3]抗小鼠IgE epsilon chain (Biotin)预吸附二抗
宿主	Rat
靶标种属	Mouse
特异性	ab11580 recognises the murine epsilon immunoglobulin heavy chain, and does not cross react with other murine immunoglobulin classes or subclasses.
经测试应用	适用于: ELISA
最小交叉反应	Rat more details
免疫原	Purified mouse IgE
偶联物	Biotin

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium azide Constituents: PBS, 50% Glycerol
纯度	Protein G purified
纯化说明	Purified IgG prepared by affinity chromatography.
克隆	单克隆
克隆编号	LO-ME-3
骨髓瘤	IR983F
同种型	IgG1
常规说明	This product should be stored undiluted. Storage in frost free freezers is not recommended. Should this product contain a precipitate we recommend microcentrifugation before use.

The Abpromise guarantee **Abpromise™**承诺保证使用ab11580于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use a concentration of 0.5 µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors