

# Rat monoclonal [LO-MM-3] Anti-Mouse IgM heavy chain ab11658

## 概述

产品名称	大鼠单克隆抗体[LO-MM-3]抗小鼠IgM heavy chain
宿主	Rat
靶标种属	Mouse
特异性	ab11658 recognises murine IgM (u heavy chain) and does not cross react with other murine immunoglobulin classes or subclasses. Avidity = $2 \times 10^8 \text{M}^{-1}$
经测试应用	适用于: ELISA
免疫原	Polyclonal mouse IgM.

## 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Purified IgM
纯化说明	Purified IgM prepared by affinity chromatography on mouse IgM/Sepharose 4B.
克隆	单克隆
克隆编号	LO-MM-3
骨髓瘤	IR983F
同种型	IgM
常规说明	This product should be stored undiluted. Storage in frost free freezers is not recommended. Should this product contain a precipitate we recommend microcentrifugation before use. This antibody recognises murine IgM ( $\mu$ heavy chain) and does not cross react with other murine immunoglobulin classes or subclasses. Avidity = $2 \times 10^8 \text{M}^{-1}$

## 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab11658于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		

#### 应用说明

ELISA: Use at a concentration of 5 µg/ml (coating).

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors