

Rat monoclonal Anti-Mouse IgG2a heavy chain [LO-MG2a-9] ab11584

概述

产品名称	大鼠单克隆抗体抗小鼠IgG2a heavy chain [LO-MG2a-9]
宿主	Rat
靶标种属	Mouse
特异性	ab11584 recognises mouse gamma 2a heavy chain. Avidity = $2.5 \times 10^9 M^{-1}$.
经测试应用	适用于: ELISA
免疫原	Purified IgG2a from BALB/c mice.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Protein G purified
纯化说明	Purified IgG prepared by affinity chromatography.
克隆	单克隆
克隆编号	LO-MG2a-9
骨髓瘤	IR983F
同种型	IgG1
常规说明	This product should be stored undiluted. Storage in frost free freezers is not recommended. Should this product contain a precipitate we recommend microcentrifugation before use.

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab11584于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use a concentration of 5 µg/ml. (coating antibody).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors