

Rat monoclonal [SB77e] Anti-Mouse IgG1 (Biotin), pre-adsorbed ab99604

2 References 1 图像


概述	
产品名称	大鼠单克隆抗体[SB77e]抗小鼠IgG1 (Biotin),预吸附二抗
宿主	Rat
靶标种属	Mouse
特异性	ab99604 reacts with Mouse IgG1 as demonstrated by ELISA. Minimal reaction to Human, Goat and Rabbit IgG. Each batch of Rat monoclonal [SB77e] Anti-Mouse IgG1 (Biotin), pre-adsorbed is tested to show minimal cross reactivity with IgM, IgA, and other subsets of IgG.
经测试应用	适用于: ELISA, Flow Cyt
最小交叉反应	Goat, Human, Rabbit
偶联物	Biotin
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. Store In the Dark.
存储溶液	pH: 7.4 Preservative: 0.1% Sodium azide Constituent: PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	SB77e
同种型	IgG2b
应用	

The Abpromise guarantee      **Abpromise™**承诺保证使用ab99604于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		1/5000 - 1/20000.
Flow Cyt		Use at an assay dependent dilution.

## 图片



Top: (Negative control) ab99604 staining for Mouse IgG1 in Human blood at  $1\mu\text{g}/10^6$  cells.

Bottom: ab99604 staining in Human blood incubated with mouse anti-Human CD3 (IgG1 isotype) and Streptavidin-FITC.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors