



Product datasheet

Rat monoclonal [SB77e] Anti-Mouse IgG1 H&L (HRP) preadsorbed ab99603

5 References

概述	
产品名称	大鼠单克隆抗体[SB77e]抗小鼠IgG1 H&L (HRP)预吸附二抗
宿主	Rat
靶标种属	Mouse
特异性	ab99603 has been tested by ELISA against a panel of purified Mouse Immunoglobulins including IgM, IgG1, IgG2a, IgG2b, IgG3, and IgA. This product is positive for IgG1 and has minimal reactivity with the other Mouse Immunoglobulins.
经测试应用	适用于: ELISA
最小交叉反应	Goat, Human, Rabbit
偶联物	HRP

[more details](#)

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.40 Constituents: 50% Glycerol (glycerin, glycerine), 50% PBS
纯度	IgG fraction
克隆	单克隆
克隆编号	SB77e
同种型	IgG2b

应用	
The Abpromise guarantee	Abpromise™ 承诺保证使用ab99603于以下的经测试应用
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。	

应用	Ab评论	说明
ELISA		1/4000 - 1/8000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors