

Sheep Anti-Rabbit IgG H&L (HRP) ab6795

[24 References](#) [2 图像](#)

概述

产品名称	羊抗兔IgG H&L (HRP)
宿主	Sheep
靶标种属	Rabbit
经测试应用	适用于: WB, ELISA, Immunomicroscopy, Dot blot, ICC/IF, IHC-P, IHC-Fr
免疫原	Rabbit IgG, whole molecule
偶联物	HRP

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.01% Gentamicin sulphate Constituents: 0.878% Sodium chloride, 1% BSA, 0.424% Potassium phosphate
纯度	Affinity purified
纯化说明	Horseradish Peroxidase (HRP).
共轭说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose beads
克隆	多克隆

应用


The Abpromise guarantee **Abpromise™** 承诺保证使用ab6795于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use at an assay dependent dilution.
ELISA		Use at an assay dependent dilution.


应用	Ab评论	说明
Immunomicroscopy		Use at an assay dependent concentration.
Dot blot		Use at an assay dependent concentration.
ICC/IF		1/500 - 1/2500.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.

图片



IHC - Wholemount - Sheep Anti-Rabbit IgG H&L (HRP) (ab6795)

This image is courtesy of an Abreview submitted by Michal Szlis



Copyright 2006 Abcam plc

Western blot - Sheep Anti-Rabbit IgG H&L (HRP) (ab6795)

ab6795 was used at dilution 1/2000 with the primary antibody **ab8895** in WB. See the review on **ab8895**.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors