

THP-1 nuclear extract lysate ab14900

概述

产品名称	THP-1核extract裂解物
常规说明	Cell line: THP1 (Human acute monocytic leukemia). Extracts have been quality control tested by Western blot and the Electrophoretic Mobility Shift Assay (EMSA).

性能

Mycoplasma free	Yes
形式	Liquid
存放说明	Shipped on Dry Ice. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.90 Constituents: 0.75% Potassium chloride, 0.009% PMSF, 0.01% Magnesium chloride, 0.008% DTT, 0.48% HEPES, 20% Glycerol
背景	Properties: differentiation; lysozyme synthesis; phagocytosis; Fc receptor; IL-1 prod.; C3b receptor. The cell line can differentiate into macrophage-like cells (phagocytic) and can be used for induction of differentiation studies. When stimulated with PMA cells have been reported to differentiate into a monocytic pathway and release arachidonic acid and prostanoids.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors